

TECHNICAL SHEET No. 16

Virus Detection: *Potato leaf roll virus* (PLRV)

Method: Immunocapture RT-PCR

General

Virus detected: PLRV from potato leaf and potato tubers.

General method was immunocapture RT-PCR

Developed by

Name of researcher: Fouad Akad and Hanokh Czosnek, The Hebrew University

Address (Email): akad@agri.huji.ac.il and Czosnek@agri.huji.ac.il

Date: December 26, 2005

Goals

To develop a sensitive method for PLRV detection based on immunocapture RT-PCR.

Methods

Immunocapture IC

1. Add 200 μ l of anti-PLRV antibody (diluted 1:1000 in coating buffer) to each ELISA plate well or to PCR tube and incubate at 37°C for 3 h. We used a commercial antibody (Bioreba).
2. Empty the ELISA well/ PCR tube and wash 3 times with TBST (10 mM Tris-HCl, pH 8.0, 150 mM NaCl, 0.05% Tween-20).
3. Homogenize 25 mg of plant leaf or tuber with 5 ml of suitable extraction buffer.
4. Add 200 μ l of the homogenate to the coated ELISA well/ PCR tube and incubate for 18 h at 4 °C.
5. Empty the ELISA well/ PCR tube and wash 3 times with TBST.
6. Dry the ELISA well/ PCR tube, add 5 to 10 μ l ddH₂O and heat at 70 °C for 15 min.

RT-PCR

cRNA synthesis

1. To 5 μ l IC add 8 μ l ddH₂O and 1 μ l primer C582 (100 pmoles); incubate at 70 °C for 10 min, followed by 30 min in an ice bath.
2. Add 1 μ l of each dNTPs (25mM each), 4 μ l reverse transcriptase buffer 5X, 1 μ l AMV reverse transcriptase (Promega); incubate at 42 °C for 1h.
3. Heat for 10 min at 90 °C; adjust the volume to 50 μ l with ddH₂O

PCR

The primers used were as follows (derived from sequence GenBank accession numbers D13954 and D13953):

C582: 5'TATCGTCCATGGGTACGGTTCGTGGT3'

R732: 5'TCTAGATCTTTGGGTTTTGCAAAGC3'

Following PCR, a 640 bp band is obtained.

1. The PCR reaction contains 5 μ l from the reverse transcriptase reaction, 0.25 μ l 25mM dNTPs, 1 μ l each primers (100 pmoles each), 2.5 μ l Taq 10 x buffer and 1 unit *Taq* polymerase; add ddH₂O to a final volume of 25 μ l.

2. Cycling: One cycle: 95 °C for 3 min, 50 °C for 2 min, 72 °C for 2 min. Thirty cycles: 95 °C for 1min, 55 °C for 1min, 72 °C for 1 min. One additional cycle: 72 °C for 10 min.
3. Analysis of reaction products: Subject the reaction products to 1% agarose gel electrophoresis.

Results

The following figure shows the detection of PLRV by IC-RT-PCR and by RT-PCR. For RT-PCR (lanes 1 – 3), RNA was from a *Datura stramonium* infected plant. For IC-RT-PCR (lanes 1' - 3'), samples were from a *Datura stramonium* infected plant. M: molecular weight marker; P: cloned PLRV DNA; 0: no template.

