

TECHNICAL SHEET No. 28

Virus detection: *Prunus necrotic ringspot virus*

Method: ELISA

General

Virus detected: *Prunus necrotic ringspot Ilarvirus* (PNRSV)

General method: ELISA

Developed by

Name of researcher: Hadia Kawtharani and Yusuf Abou-Jawdah, American University of Beirut,
Lebanon

Address (email): haddoush21@hotmail.com and abujawyf@aub.edu.lb

Date: November 23, 2002

Goals

Identification of PNRSV strains and production of PNRSV ELISA Kits.

Introduction

Prunus necrotic ringspot virus

Phosphate Buffer Saline (PBS)	Washing Buffer pH 7.4	Coating Buffer pH 9.6
Distilled water 1.0	PBS 1.0 l	Distilled water 1.0 l
NaCl 8.0 g	TWEEN 20 0.5 ml	Na ₂ CO ₃ 1.6 g
Na ₂ HPO ₄ .12H ₂ O 2.9 g		NaHCO ₃ 3.0 g
KH ₂ PO ₄ 0.2 g		
KCl 0.2 g		

Conjugate Buffer pH 7.1-7.3	Substrate Buffer pH 9.6 adjusted with concentrated HCl
PBS TWEEN 100 ml	Distilled water 80.0 ml
PVP 2 g	Diethanolamine 9.7 ml
Ovalbumine 0.2 g	

Coating of ELISA plates

1. Dilute antibodies in coating buffer as recommended by the supplier.
2. Add 100 µl to each well, cover the plate.
3. Incubate for 3-4 hr at 37°C.
4. Wash three times with PBS-TWEEN.

Antigen extraction and binding

1. Leaf tissues (1:10 w/v) were extracted in 0.1 M phosphate buffer Saline-Tween-Polyvinylpyrrolidone (PBST-PVP), pH 7.4 containing 2% PVP, 0.05% Tween-20, 0.15 M NaCl, and 4 mM KCl.
2. Add 100 µl extracted sample/well; two wells were used per sample.
3. Incubate overnight at 4°C.
4. Wash three times with PBS-TWEEN.

Conjugate

1. Dilute conjugated antibodies in conjugate buffer as recommended and add 100 µl to each well.
2. Incubate for 3-4 h at 37°C.
3. Wash three times with PBS-TWEEN.

Color reaction

1. Dissolve 10 mg pNPP in 10 ml substrate buffer just before use.
2. Incubate for 1 h.
3. After incubation, the reaction was detected calorimetrically at A₄₀₅ nm using an ELISA reader.
4. The test was considered positive when the mean absorbance value of a sample was over twice that of healthy controls.